


GOAT ANTI-RNF35 / TRIM40 ANTIBODY

SKU: EB05890



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names TRIM40|tripartite motif-containing 40|RNF35|ring finger RNF35|DAAP-200B17.3|ring finger protein 40

Accession ID NP_619645.1

Blocking Peptide EBP05890

Immunogen Peptide with sequence C-GVSELLLPQKL, from the C Terminus of the protein sequence according to NP_619645.1.

Peptide Sequence C-GVSELLLPQKL

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species Human

Reactive Species Human

Human Gene ID 135644

Product Grade https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection Limit Antibody detection limit dilution 1:32000.

Western Blot Approx 26kDa band observed in Human Placenta lysates (calculated MW of 26.1kDa according to NP_619645.1). Recommended concentration: 0.5-2µg/ml.

Application Type Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa