

GOAT ANTI-RNF7 / SAG / ROC2 ANTIBODY

SKU: EB05842



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	zinc RING finger protein SAG sensitive to apoptosis, zinc RING finger protein SAG, regulator of cullins 2 regulator of cullins 2 sensitive to apoptosis gene ring finger protein 7 CKBBP1 ROC2 SAG RNF7
Accession ID	NP_055060.1; NP_001188299.1
Blocking Peptide	EBP05842
Immunogen	Peptide with sequence CQQDWVVQRIGK, from the C Terminus of the protein sequence according to NP_055060.1; NP_001188299.1.
Product Comments	This antibody is expected to recognise isoform 1 (NP_055060.1) and isoform 4 (NP_001188299.1).
Peptide Sequence	CQQDWVVQRIGK
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	9616
Mouse Gene ID	19823
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	No signal obtained yet but low background observed in Human Heart and Human Testis lysates at upto 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)