

GOAT ANTI-ROBO3 / RIG1 (INTERNAL) ANTIBODY

SKU: EB07773

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	roundabout, axon guidance receptor, homolog 3 retinoblastoma inhibiting gene
Names	1 RBIG1 HGPPS FLJ21044 roundabout, axon guidance receptor, homolog 3 (Drosophila) RIG1 ROBO3
Accession ID	NP_071765.2
Blocking Peptide	EBP07773
Immunogen	Peptide with sequence C-RRQRKELSHYT, from the internal region of the protein sequence according to NP_071765.2.
Peptide Sequence	C-RRQRKELSHYT
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Pig, Cow
Reactive Species	Human
Human Gene ID	64221
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 150kDa band observed in Human Brain (Cerebellum and Hippocampus) lysates (calculated MW of 148kDa according to NP_071765.2). Recommended concentration: 0.5-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

