


GOAT ANTI-RP11 / PRPF31 ANTIBODY

SKU: EB05363



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	PRP31 NY-BR-99 PRPF31 PRP31 pre-mRNA processing factor 31 homolog (yeast) RP11 DKFZp566j153 pre-mRNA processing factor 31 homolog (yeast) Protein 61K PRP31 pre-mRNA processing factor 31 homolog (S. cerevisiae) pre-mRNA processing factor 31 homolog
Accession ID	NP_056444.3
Blocking Peptide	EBP05363
Immunogen	Peptide with sequence C-LKVKGEKSGLMST, from the C Terminus of the protein sequence according to NP_056444.3.
Peptide Sequence	C-LKVKGEKSGLMST
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog
Reactive Species	Human
Human Gene ID	26121
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Approx 60kDa band observed in lysates of cell line A431 (calculated MW of 55.4kDa according to NP_056444.3). Recommended concentration: 0.3-1.0µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

