

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RPGRIP1L ANTIBODY

SKU: EB07988



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias MKS5|KIAA1005|JBTS7|DKFZp686C0668|CORS3|RPGRIP1-like

Names

Accession

NP_056087.2 ID

Blocking

EBP07988 **Peptide**

Peptide with sequence C-DKENNKAKRDILK, from the internal region of the protein sequence according to **Immunogen**

NP 056087.2.

Peptide Sequence

C-DKENNKAKRDILK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Dog **Species**

Reactive

Human **Species**

Human

23322 **Gene ID**

Product https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Western Approx 150kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 151kDa according to

Blot NP_056087.2). Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB Type

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa