



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RPS19 ANTIBODY

SKU: EB07023



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

40S ribosomal protein S19|DBA|HGNC:10402|ribosomal protein S19|RPS19

Names

Accession

NP_001013.1

Blocking

EBP07023

Peptide

Peptide with sequence C-KMVEKDQDGGRK, from the internal region of the protein sequence according to **Immunogen**

NP 001013.1.

Peptide

C-KMVEKDQDGGRK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat, Dog, Cow **Species**

Reactive

Human **Species**

Human

6223 Gene ID

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. a mazon aws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Brain (Cortex), Liver and Tonsil. Recommended concentration: 5µg/ml.

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Blot

Approx 16kDa band observed in lysate of cell line K562 (calculated MW of 16.1kDa according to NP 001013.1). Western

Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.

Application

Pep-ELISA, WB, IHC **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa


50kDa


37kDa

25kDa

20kDa

15kDa









Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

