








GOAT ANTI-RSK2 ANTIBODY

SKU: EB06985



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / ribosomal protein S6 kinase, 90kD, polypeptide 3|mental retardation, X-linked 19|insulin-stimulated protein kinase 1|RSK|CLS|RP11-393H10.3|pp90RSK2|p90-RSK2|S6K-alpha3|RSK2|MRX19|MAPKAPK1B|ISPK-1|HU-3|HGNC:10432|ribosomal

protein S6 kinase, 90kDa, polypeptide 3|RPS6KA3

Names Accession

Alias

ID

NP_004577.1

Immunogen Peptide with sequence C-DDESQAMQTVGVH, from the internal region of the protein sequence according to NP 004577.1.

Peptide Sequence

C-DDESQAMQTVGVH

Method

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human, Mouse, Rat

Species Reactive

Human

Species Human

Gene ID

6197





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Product

 $https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png$

Grade

 $\textbf{IHC Results} \ \ \text{Paraffin embedded Human Skin. Recommended concentration: } 5 \mu g/ml.$

ELISA

Detection Antibody detection limit dilution 1:64000.

Limit

Application Pep-ELISA, IHC Type

GALLERY IMAGES

