



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-RSPO3 ANTIBODY

SKU: EB10509



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Unit Size

100 µg

Synonym / **Alias**

RSPO3|THSD2|thrombospondin-like|thrombospondin, type I, domain containing 2|R-spondin 3 homolog

(Xenopus laevis)|R-spondin 3|PWTSR|FLJ14440|CRISTIN1

Names

Accession NP_116173.2

Blocking

EBP10509

Peptide

Immunogen

Peptide with sequence QHPSAKGNLCPPT, from the internal region of the protein sequence according to

NP 116173.2.

Peptide

Sequence

QHPSAKGNLCPPT

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Cow

Reactive

Species

Human Gene ID

84870

Human

Mouse

72780 **Gene ID**

Rat Gene ID 498997

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

ELISA

Antibody detection limit dilution 1:8000. Detection

Limit

Western Approx 28kDa band observed in Human Skeletal Muscle lysates (calculated MW of 30.9kDa according to

Blot NP_116173.2). Recommended concentration: 0.3-1μg/ml. Primary incubation was 1 hour.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa