

GOAT ANTI-RSPO3 ANTIBODY

SKU: EB10509

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	RSPO3 THSD2 thrombospondin-like thrombospondin, type I, domain containing 2 R-spondin 3 homolog (Xenopus laevis) R-spondin 3 PWTSR FLJ14440 CRISTIN1
Accession ID	NP_116173.2
Blocking Peptide	EBP10509
Immunogen	Peptide with sequence QHPSAKGNLCPPT, from the internal region of the protein sequence according to NP_116173.2.
Peptide Sequence	QHPSAKGNLCPPT
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Cow
Reactive Species	Human
Human Gene ID	84870
Mouse Gene ID	72780
Rat Gene ID	498997
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 28kDa band observed in Human Skeletal Muscle lysates (calculated MW of 30.9kDa according to NP_116173.2). Recommended concentration: 0.3-1µg/ml. Primary incubation was 1 hour.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

