Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SAFB2 ANTIBODY

SKU: EB09246

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias KIAA0138|scaffold attachment factor B2|SAFB2

Names

Accession NP 055464.1

Blocking

EBP09246 **Peptide**

Peptide with sequence C-GEHNQRLEEHQARA, from the internal region of the protein sequence according to Immunogen

NP 055464.1.

Product Comments

This antibody is not expected to cross-react with SAFB1.

Peptide

C-GEHNQRLEEHQARA Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Human

9667 **Gene ID**

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:32000.

Limit

Application Pep-ELISA **Type**

> Goat Anti-SAFB2 Antibody https://everestbiotech.com/products/goat-anti-safb2-antibody/