

GOAT ANTI-SAFB2 ANTIBODY

SKU: EB09246

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	KIAA0138 scaffold attachment factor B2 SAFB2
Names	
Accession ID	NP_055464.1
Blocking Peptide	EBP09246
Immunogen	Peptide with sequence C-GEHNQRLEE HQARA, from the internal region of the protein sequence according to NP_055464.1.
Product Comments	This antibody is not expected to cross-react with SAFB1.
Peptide Sequence	C-GEHNQRLEE HQARA
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	9667
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)