

GOAT ANTI-SAK (STK18) ANTIBODY

SKU: EB05225



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Snk Akin Kinase Sak SAK serine/threonine kinase 18 STK18
Accession ID	NP_055079
Blocking Peptide	EBP05225
Immunogen	Peptide with sequence C-SILLMFSNPTPNFH, from the C Terminus of the protein sequence according to NP_055079.
Peptide Sequence	C-SILLMFSNPTPNFH
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	10733
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Human Brain and Kidney extracts at upto 1µg/ml.