


Telephone: (650) 697-3600

## **GOAT ANTI-SCARB1 / SR-BI (AA238-250) ANTIBODY**

**SKU:** EB12301



## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / scavenger receptor class B type III|scavenger receptor class B member 1|CD36 antigen (collagen type I

receptor, thrombospondin receptor)-like 1|CD36 and LIMPII analogous 1|SRB1|SR-Alias **Names** BI|HDLQTL6|CLA1|CLA-1|CD36L1|scavenger receptor class B, member 1|SCARB1

**Accession** 

ID

NP\_005496.4; NP\_001076428.1

**Blocking** 

**Peptide** 

EBP12301

**Immunogen** 

Peptide with sequence C-NGLSKVDFWHSDQ, from the internal region of the protein sequence according to

NP\_005496.4; NP\_001076428.1.

**Product Comments** 

This antibody is expected to recognize both reported isoforms (NP 005496.4; NP 001076428.1).

**Peptide** 

Sequence

C-NGLSKVDFWHSDQ

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

Shipping Instructions

Refrigerated

**Predicted** 

**Species** 

Human

Human

Reactive

**Species** 

Human **Gene ID** 

949

**Product** 

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite\_medium.png$ Grade

IHC Results Paraffin embedded Human Testis. Recommended concentration: 5µg/ml.

**ELISA** 

Detection Antibody detection limit dilution 1:128000.

Limit

Western Blot

Approx 85kDa band observed in Human Adrenal gland lysates (calculated MW of 57.0kDa according to

NP\_005496.4). Recommended concentration: 0.01-0.03µg/ml.

**Type** 


**Application** Pep-ELISA, WB, IHC

## **GALLERY IMAGES**









250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa