

GOAT ANTI-SCRIB / SCRIBBLE ANTIBODY

SKU: EB09387



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage	Aliquot and store at -20°C. Minimize freezing and thawing.
Instructions	
Synonym /	
Alias	scribble Vartul SCRIB1 SCRB1 CRIB1 scribbled homolog (Drosophila) SCRIB
Names	
Usage	Immunofluorescence: Strong expression of the protein seen in the cell junctions of MCF7
Summary	cells. Recommended concentration: 10µg/ml.
Accession	
ID	NP_874365.3; NP_056171.3
Blocking	
Peptide	EBP09387
Immunogen	Peptide with sequence C-PEGPGKEKELPGQ, from the internal region of the protein sequence according to NP_874365.3; NP_056171.3.
Product	
Comments	This antibody is expected to recognize both reported isoforms (NP_874365.3; NP_056171.3).
Peptide	
Sequence	C-PEGPGKEKELPGQ
Purification	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography
Method	using the immunizing peptide.
Shipping	
Instructions	Refrigerated
Predicted	
Species	Human
Reactive	
Species	Human
Human	
Gene ID	23513
Product	
Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection	Antibody detection limit dilution 1:1000.
Limit	
Western	Approx 200-250kDa band observed in lysates of cell lines HeLa, HepG2 and U2OS (calculated MW of 178kDa
Blot	according to NP_874365.3). This molecular weight is routinely observed by other sources. Recommended
	concentration: 1-3µg/ml. Primary incubation 1 hour at room temperature.
Application	
Type	Pep-ELISA, WB, IF

GALLERY IMAGES

