


GOAT ANTI-SCRIBBLE / VARTUL ANTIBODY

SKU: EB06251



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names scribble|KIAA0147|Vartul|SCRIB1|SCRB1|CRIB1|Circletail homolog (Mouse)|scribbled homolog (Drosophila)|SCRIB

Accession ID NP_874365.2; NP_056171.2

Blocking Peptide EBP06251

Immunogen Peptide with sequence C-VRPGRRGLGPVPS, from the C Terminus of the protein sequence according to NP_874365.2; NP_056171.2.

Peptide Sequence C-VRPGRRGLGPVPS

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human, Mouse, Rat
Reactive Species	Human
Human Gene ID	23513
Mouse Gene ID	105782
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Tonsil. Recommended concentration: 6-8µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:8000.
Application Type	Pep-ELISA, IHC

GALLERY IMAGES

