
Telephone: (650) 697-3600





GOAT ANTI-SCRIBBLE / VARTUL ANTIBODY

SKU: EB06251



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

scribble|KIAA0147|Vartul|SCRIB1|SCRB1|CRIB1|Circletail homolog (Mouse)|scribbled homolog (Drosophila)|SCRIB

Alias Names

Accession ID

NP_874365.2; NP_056171.2

Blocking Peptide

EBP06251

Immunogen

Peptide with sequence C-VRPGRRGLGPVPS, from the C Terminus of the protein sequence according to NP_874365.2; NP_056171.2.

Peptide

C-VRPGRRGLGPVPS

Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human, Mouse, Rat

Species

Reactive **Species**

Human

Human

23513

Gene ID Mouse

Gene ID

105782

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Tonsil. Recommended concentration: 6-8µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:8000.

Limit

Type

Application Pep-ELISA, IHC

GALLERY IMAGES

