


GOAT ANTI-SDCCAG8 (AA469-483) ANTIBODY

SKU: EB11472



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / Alias Names antigen NY-CO-8|BBS16|CCCAP|centrosomal colon cancer autoantigen protein|hCCCAP|NPHP10|NY-CO-8|OTTHUMP00000037913|SDCCAG8|serologically defined colon cancer antigen 8|SLSN7

Accession ID NP_006633.1

Blocking Peptide EBP11472

Immunogen Peptide with sequence C-DEAEKEHREFRAKTN, from the internal region of the protein sequence according to NP_006633.1.

Peptide Sequence C-DEAEKEHREFRAKTN

Purification Method Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.

Shipping Instructions Refrigerated

Predicted Species	Human
Reactive Species	Human
Human Gene ID	10806
Mouse Gene ID	76816
Rat Gene ID	305002
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
IHC Results	Paraffin embedded Human Brain (Cortex). Recommended concentration: 5µg/ml.
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Preliminary experiments in Human Prostate, Tonsil and Thymus lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Application Type	Pep-ELISA, IHC

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

