


Telephone: (650) 697-3600

GOAT ANTI-SDHAF1 ANTIBODY

SKU: EB09604



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

SDHAF1, succinate dehydrogenase complex assembly factor 1

Names

Accession

NP_001036096.1

Blocking

ID

EBP09604

Peptide

Peptide with sequence C-HDSTGAPETRPDGR, from the C Terminus of the protein sequence according to Immunogen

NP_001036096.1.

Peptide Sequence

C-HDSTGAPETRPDGR

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping Instructions

Refrigerated





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human, Mouse, Dog

Species Reactive

Human

Species

Human **Gene ID**

644096

Mouse **Gene ID**

68332

Product

 $https://prod-vector-labs-pimcore-assets.s 3. us-east-1. amazonaws.com/assets/products/image/elite_medium.png$ Grade

IHC Results Paraffin embedded Human Placenta. Recommended concentration: 5µg/ml.

ELISA

Detection Antibody detection limit dilution 1:16000.

Limit

Approx 19kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 12.8kDa according to Western

Blot NP_001036096.1). Recommended concentration: $0.3-1\mu g/ml$. Primary incubation was 1 hour.


Application Pep-ELISA, WB, IHC **Type**

GALLERY IMAGES



Telephone: (650) 697-3600





250kDa 150kDa 100kDa 75kDa 50kDa 37kDa 25kDa 20kDa 15kDa