



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SEC61A1 ANTIBODY

SKU: EB05001







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias Names Sec61 alpha-1|Sec61 alpha 1 subunit|Sec61 alpha 1 subunit (S. cerevisiae)|sec61 homolog|protein transport

protein SEC61 alpha subunit isoform 1|SEC61A|SEC61|HSEC61|SEC61A1

Accession

ID

NP_037468.1

Blocking

Peptide

EBP05001

Peptide with sequence C-KEQSEVGSMGALLF, from the C Terminus of the protein sequence according to **Immunogen**

NP 037468.1.

Peptide

Sequence

C-KEQSEVGSMGALLF

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Dog, Cow, Pig

Reactive

Species

Human

Human

29927 Gene ID

Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

ELISA

Detection

Antibody detection limit dilution 1:16000.

Limit

Western **Blot**

Approx 48kDa band observed in Human Brain (Cerebelum) lysates (calculated MW of 52.3kDa according to

NP_037468.1). Recommended concentration: 0.3-1µg/ml.

Application

Type

Pep-ELISA, WB

GALLERY IMAGES







250kDa 150kDa 100kDa 75kDa

50kDa

37kDa

25kDa

20kDa

15kDa