

GOAT ANTI-SECRETOGRANIN II (AA203-217) ANTIBODY

SKU: EB12272



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	secretoneurin secretogranin-2 chromogranin-C SgII SN EM66 CHGC secretogranin II SCG2
Accession ID	NP_003460.2
Blocking Peptide	EBP12272
Immunogen	Peptide with sequence C-QELGKLTGPNNQKRE, from the internal region of the protein sequence according to NP_003460.2.
Peptide Sequence	C-QELGKLTGPNNQKRE
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Dog, Pig
Reactive Species	Human
Human Gene ID	7857
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 70kDa band observed in Human Brain (Cerebellum) lysates (calculated MW of 70.9kDa according to NP_003460.2). Recommended concentration: 0.3-1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

