

GOAT ANTI-SELECTIN E ANTIBODY

SKU: EB12043



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias	leukocyte-endothelial cell adhesion molecule 2 endothelial leukocyte adhesion molecule 1 endothelial adhesion molecule 1 ELAM-1 E-selectin CD62 antigen-like family member
Names	E LECAM2 ESEL ELAM1 ELAM CD62E selectin E SELE
Accession ID	NP_000441.2
Blocking Peptide	EBP12043
Immunogen	Peptide with sequence SPEHGSLVCSHP, from the internal region of the protein sequence according to NP_000441.2.
Peptide Sequence	SPEHGSLVCSHP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	6401
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Preliminary experiments in Human Heart, Lung, Umbilical Cord lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml).
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)