

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SELECTIN E ANTIBODY

SKU: EB12043

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / leukocyte-endothelial celleukocyte endothelial cell adhesion molecule 2|endothelial leukocyte adhesion molecule

Alias 1|endothelial adhesion molecule 1|ELAM-1|E-selectin|CD62 antigen-like family member

Names E|LECAM2|ESEL|ELAM1|ELAM|CD62E|selectin E|SELE

Accession

NP 000441.2 ID

Blocking

EBP12043 **Peptide**

Peptide with sequence SPEHGSLVCSHP, from the internal region of the protein sequence according to Immunogen

NP 000441.2.

Peptide

SPEHGSLVCSHP Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Species

Human

Human

6401 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection Antibody detection limit dilution 1:128000.

Limit

Blot

Western Preliminary experiments in Human Heart, Lung, Umbillical Cord lysates gave no specific signal but low background

(at antibody concentration up to 1µg/ml).

Application

Pep-ELISA **Type**