

GOAT ANTI-SENP1 ANTIBODY

SKU: EB05334



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	sentrin/SUMO-specific protease SENP1
Usage Summary	raju
Accession ID	NM_014554
Immunogen	Peptide with sequence DDIADRMAMDAGEV-C, from the N Terminus of the protein sequence according to NM_014554.
Peptide Sequence	DDIADRMAMDAGEV-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	29843
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Human Brain, Mouse Liver and Daudi cell lysates at upto 2µg/ml.