

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SENP1 ANTIBODY**

**SKU:** EB05334

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias sentrin/SUMO-specific protease|SENP1

**Names** 

**Usage** raju **Summary** 

Accession

NM\_014554 ID

Immunogen NM\_014554. Peptide with sequence DDIADRMRMDAGEV-C, from the N Terminus of the protein sequence according to

**Peptide** 

DDIADRMRMDAGEV-C Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

**Species** 

Human

Human **Gene ID** 

29843

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Western

Western blot: No signal obtained yet but low background observed in Human Brain, Mouse Liver and Daudi cell

Blot lysates at upto 2µg/ml.