

GOAT ANTI-SENP7 (N TERMINUS) ANTIBODY

SKU: EB06090



SPECIFICATIONS

Unit Size	100 µg
Synonym /	
Alias	KIAA1707 sentrin/SUMO-specific protease SENP7
Names	
Accession ID	NP_065705
Blocking Peptide	EBP06090
Immunogen	Peptide with sequence DKRKLGRPPSSSE-C, from the N Terminus of the protein sequence according to NP_065705.
Peptide Sequence	DKRKLGRPPSSSE-C
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	57337
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:16000.
Western Blot	Western Blot: Preliminary experiments gave no signal but low background in Human Testis and Kidney extracts at up to 1µg/ml. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA