


Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SEPT6 ANTIBODY





SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

septin 2|SEPT2|SEP2|MGC20339|MGC16619|KIAA0128|septin 6|SEPT6 **Alias**

Names

Immunofluorescence: Strong expression of the protein seen in A431 and U2OS cells.

Usage Recommended concentration: 10µg/ml. Flow Cytometry: Flow cytometric analysis of A431 Summary

cells. Recommended concentration: 10ug/ml.

Accession

NP_665798.1; NP_055944.2; NP_665801.1

Blocking

ID

EBP08932 **Peptide**





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Peptide with sequence C-DEVNAFKQRKTA, from the internal region of the protein sequence according to **Immunogen**

NP_665798.1; NP_055944.2; NP_665801.1.

Product This antibody is expected to recognize the reported isoforms NP 665798.1, NP 055944.2 and NP 665801.1.

Reported variants NP_665799.1 and NP_665798.1 represent identical protein. Comments

Peptide C-DEVNAFKQRKTA Sequence

e erest

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human, Cow

Species

Reactive Human **Species**

Human

23157 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite plus medium.png Grade

In paraffin embedded Human Tonsil shows specific staining in dividing cells just outside a germinal centre. **IHC Results**

Recommended concentration: 3-5µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:4000.

Limit Western

Approx 50kDa band observed in Human Testis, Duodenum, Kidney and Tonsil lysates, and in lysates of cell lines Daudi, Jurkat and MOLT-4 (calculated MW of 49.7kDa according to NP_055944.2). Recommended concentration:

0.1-1µg/ml. Primary incubation 1 hour at room temperature.


Application


Туре


Blot

Pep-ELISA, WB, IHC, IF, FC

GALLERY IMAGES











Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

