Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SERPINA10 / ZPI ANTIBODY**

**SKU:** EB08624

## **SPECIFICATIONS**

**Unit Size** 100 µg

Synonym / serine (or cysteine) proteinase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 10|protein Z-Alias dependent protease inhibitorprecursor|protein Z-dependent protease inhibitor|PZI|serpin peptidase inhibitor,

**Names** clade A (alpha-1 antiproteinase, antitrypsin), member 10|ZPI|SERPINA10

Accession

NP 001094077.1; NP 057270.1 ID

**Blocking** EBP08624 **Peptide** 

Peptide with sequence QSPETPAPQNQTSR-C, from the N Terminus (near) of the protein sequence according to **Immunogen** 

NP\_001094077.1; NP\_057270.1.

This antibody is expected to recognize the isoforms (NP 001094077.1; NP 057270.1) that are shwoing 100% **Product** 

Comments sequence identity.

**Peptide** 

QSPETPAPQNQTSR-C Sequence

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human

**Species** 

Human 51156

**Gene ID** 

**Product** 

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:16000.

Limit

Application

Pep-ELISA **Type**