

# GOAT ANTI-SERPINA10 / ZPI ANTIBODY

SKU: EB08624



---

## SPECIFICATIONS

<b>Unit Size</b>	100 µg
<b>Synonym / Alias</b>	serine (or cysteine) proteinase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 10 protein Z-dependent protease inhibitorprecursor protein Z-dependent protease inhibitor PZI serpin peptidase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 10 ZPI SERPINA10
<b>Accession ID</b>	NP_001094077.1; NP_057270.1
<b>Blocking Peptide</b>	EBP08624
<b>Immunogen</b>	Peptide with sequence QSPETPAPQNQTSR-C, from the N Terminus (near) of the protein sequence according to NP_001094077.1; NP_057270.1.
<b>Product Comments</b>	This antibody is expected to recognize the isoforms (NP_001094077.1; NP_057270.1) that are showing 100% sequence identity.
<b>Peptide Sequence</b>	QSPETPAPQNQTSR-C
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Reactive Species</b>	Human
<b>Human Gene ID</b>	51156
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>ELISA Detection Limit</b>	Antibody detection limit dilution 1:16000.
<b>Application Type</b>	Pep-ELISA

## DOCUMENTS

- [Data Sheet](#)

