

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SERPINA12 ANTIBODY

SKU: EB08542



Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / visceral adipose-specific SERPIN|vaspin|serine (or cysteine) proteinase inhibitor, clade A (alpha-1 antiproteinase, antitrypsin), member 12|OL-64|serpin peptidase inhibitor, clade A (alpha-1 antiproteinase,

Alias antitrypsin), member 12|SERPINA12 Names

Accession

NP_776249.1 ID

Blocking

EBP08542 **Peptide**

Peptide with sequence C-RARWKHEFDPNVTK, from the internal region of the protein sequence according to **Immunogen**

NP 776249.1.

Peptide

C-RARWKHEFDPNVTK Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

Predicted

Human **Species**

Reactive

Human **Species**

Human

145264 Gene ID

Product

Grade

ELISA

Antibody detection limit dilution 1:32000. **Detection**

Limit

Western Approx 48kDa band observed in Human Adipose lysates (calculated MW of 47.2kDa according to

Blot NP_776249.1). Recommended concentration: 0.3-1µg/ml.

Application

Pep-ELISA, WB **Type**

GALLERY IMAGES



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa