

GOAT ANTI-SERPINB11 ANTIBODY

SKU: EB08539



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SERPINB11 serpin peptidase inhibitor, clade B (ovalbumin), member 11 EPIPIN MGC163242 MGC163244 serine (or cysteine) proteinase inhibitor, clade B (ovalbumin), member 11 serpin peptidase inhibitor, clade B, member 11
Accession ID	NP_536723.2
Blocking Peptide	EBP08539
Immunogen	Peptide with sequence EFSQINQPDSNCT, from the internal region of the protein sequence according to NP_536723.2.
Peptide Sequence	EFSQINQPDSNCT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	89778
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Preliminary experiments in Human Lung, Placenta, Prostate and Tonsil lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)