

GOAT ANTI-SERPINB6 ANTIBODY

SKU: EB08577

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SERPINB6 serpin peptidase inhibitor, clade B (ovalbumin), member 6 CAP DKFZp686I04222 MGC111370 MSTP057 PI6 PTI protease inhibitor 6 (placental thrombin inhibitor) serine (or cysteine) proteinase inhibitor, clade B (ovalbumin), member 6
Accession ID	NP_004559.4
Blocking Peptide	EBP08577
Immunogen	Peptide with sequence C-NWDEQFDKENTEER, from the internal region of the protein sequence according to NP_004559.4.
Peptide Sequence	C-NWDEQFDKENTEER
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	5269
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 40kDa band observed in Human Ovary lysates (calculated MW of 42.6kDa according to NP_004559.4). Recommended concentration: 0.1-.3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

