

Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-SERPINB6 ANTIBODY**

**SKU:** EB08577



Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

250kDa 150kDa 100kDa


75kDa 50kDa

37kDa

25kDa

20kDa

15kDa





## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym / SERPINB6|serpin peptidase inhibitor, clade B (ovalbumin), member

**Alias** 6|CAP|DKFZp686I04222|MGC111370|MSTP057|PI6|PTI|protease inhibitor 6 (placental thrombin inhibitor)|serine

(or cysteine) proteinase inhibitor, clade B (ovalbumin), member 6 Names

**Accession** 

NP 004559.4 ID

**Blocking** 

EBP08577 **Peptide** 

Peptide with sequence C-NWDEQFDKENTEER, from the internal region of the protein sequence according to **Immunogen** 

NP 004559.4.

**Peptide** 

C-NWDEQFDKENTEER Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

**Shipping** 

Refrigerated Instructions

**Predicted** 

Human **Species** 

Reactive

Human **Species** 

Human

5269 Gene ID

**Product** 

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png Grade

**ELISA** 

Antibody detection limit dilution 1:128000. **Detection** 

Limit

Western Approx 40kDa band observed in Human Ovary lysates (calculated MW of 42.6kDa according to NP\_004559.4).

**Blot** Recommended concentration: 0.1-.3µg/ml.

**Application** 

Pep-ELISA, WB **Type** 

## **GALLERY IMAGES**



Telephone: (650) 697-3600



250kDa 150kDa 100kDa

75kDa 50kDa

37kDa

25kDa

20kDa

15kDa