

GOAT ANTI-SETBP1 ANTIBODY

SKU: EB05242

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	KIAA0437 Seb SEB SET binding protein 1 SETBP1
Accession ID	NP_056374
Immunogen	Peptide with sequence C-QRKSARGESEVLP, from the C Terminus of the protein sequence according to NP_056374.
Peptide Sequence	C-QRKSARGESEVLP
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Human Gene ID	26040
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Human Brain extracts at upto 2µg/ml.