






GOAT ANTI-SFRP1 ANTIBODY

SKU: EB08864



SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

Alias

secreted apoptosis-related protein 2|SARP2|FrzA|FRP1|FRP-1|FRP|secreted frizzled-related protein 1|SFRP1

Names

Accession

NP_003003.3

Blocking

EBP08864

Peptide

Peptide with sequence KNHECPTFQSVFK, from the C Terminus of the protein sequence according to

Immunogen NP 003003.3.

Peptide Sequence

KNHECPTFQSVFK

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions





Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

Predicted

Human, Mouse, Cow

Species

Reactive

Human **Species**

Human

6422

Gene ID Product

Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png

IHC Results

In paraffin embedded Human Kidney shows vesicular staining in epithelial cells of PCT. Recommended

concentration, 3-6µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:1000.

Limit

Application Pep-ELISA, IHC

Type

GALLERY IMAGES

