

GOAT ANTI-SFRP2 ANTIBODY

SKU: EB07622



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	
Names	secreted apoptosis related protein 1 SDF-5 SARP1 FRP-2 secreted frizzled-related protein 2 SFRP2
Accession ID	NP_003004.1
Blocking Peptide	EBP07622
Immunogen	Peptide with sequence C-KRWQKGQREFKR, from the C Terminus of the protein sequence according to NP_003004.1.
Peptide Sequence	C-KRWQKGQREFKR
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog
Reactive Species	Human
Human Gene ID	6423
Mouse Gene ID	20319
Rat Gene ID	310552
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA	
Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Approx 38kDa band observed in Human Duodenum lysates (calculated MW of 33.5kDa according to NP_003004.1). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

