




Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SFRP4 ANTIBODY

SKU: EB09115







250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa 20kDa

15kDa





Telephone: (650) 697-3600

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size 100 μg

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

secreted frizzled-related protein 4|MGC26498|FRPHE|FRP-4|secreted frizzled-related protein 4|SFRP4

Names

ID

Accession

NP_003005.2

Blocking

Peptide

EBP09115

Immunogen

Peptide with sequence TNSSCQCPHILPHQD, from the internal region of the protein sequence according to

NP 003005.2.

Peptide

Sequence

TNSSCQCPHILPHQD

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Instructions

Refrigerated

Predicted

Species

Human, Mouse, Rat, Cow, Pig

Reactive

Species

Human

Human

6424 Gene ID

Mouse

Gene ID

Rat Gene ID 89803

20379

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png Grade

IHC Results

In paraffin embedded Human Uterus shows granular staining in the apical cytoplasm of the endometrial gland.

Recommended concentration, 2-4µg/ml.

ELISA

Detection

Antibody detection limit dilution 1:32000.

Limit

Western Approx 40kDa band observed in Human Uterus lysates (calculated MW of 39.8kDa according to NP_003005.2).

Blot Recommended concentration: 1-3µg/ml.

Application

Type

Pep-ELISA, WB, IHC





Telephone: (650) 697-3600

GALLERY IMAGES

250kDa 150kDa 100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

