

GOAT ANTI-SH2B (683 AA ISOFORM) ANTIBODY

SKU: EB05092



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	likely ortholog of mouse SH2 domain-containing putative adapter SH2-B SH2-B gamma signaling protein SH2-B alpha signaling protein DKFZP547G1110 protein KIAA1299 DKFZP547G1110 SH2-B SH2-B homolog SH2B
Accession ID	AAF73914.1
Immunogen	Peptide with sequence MPPAPSCPSERVTV, from the C Terminus of the protein sequence according to AAF73914.1.
Product Comments	This antibody will recognise the 683 aa isoform.
Peptide Sequence	MPPAPSCPSERVTV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse
Human Gene ID	25970
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: No signal obtained yet but low background observed in Jurkat and U937 lysates at upto 2µg/ml.