

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SH2B (683 AA ISOFORM) ANTIBODY

SKU: EB05092

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

likely ortholog of mouse SH2 domain-containing putative adapter SH2-B|SH2-B gamma signaling protein|SH2-B **Alias** alpha signaling protein|DKFZP547G1110 protein|KIAA1299|DKFZP547G1110|SH2-B|SH2-B homolog|SH2B **Names**

Accession AAF73914.1

Immunogen Peptide with sequence MPPAPSCPSERVTV, from the C Terminus of the protein sequence according to AAF73914.1.

Product Comments

This antibody will recognise the 683 aa isoform.

Peptide

MPPAPSCPSERVTV Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

Shipping Instructions

Refrigerated

Predicted

Human, Mouse **Species**

Human

25970 Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png Grade

Western

Western blot: No signal obtained yet but low background observed in Jurkat and U937 lysates at upto 2µg/ml. Blot