

GOAT ANTI-SH2B (756 AA ISOFORM) ANTIBODY

SKU: EB05093

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	likely ortholog of mouse SH2 domain-containing putative adapter SH2-B SH2-B gamma signaling protein SH2-B alpha signaling protein DKFZP547G1110 protein KIAA1299 DKFZP547G1110 SH2-B SH2-B homolog SH2B
Immunogen	Peptide with sequence C-GHPRAINNQYSFV, from the C Terminus of the protein sequence according to .
Product Comments	This antibody will recognise the 756 aa isoform.
Peptide Sequence	C-GHPRAINNQYSFV
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat
Human Gene ID	25970
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
Western Blot	Western blot: Two consistent bands observed in Jurkat and U937 at approx 85kDa and approx 45kDa with 1 hour incubation of primary.