

GOAT ANTI-SH2D1A / SAP ANTIBODY

SKU: EB05103



SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias	SLAM (signaling lymphocyte activation molecule) associated protein SH2 domain protein 1A XLPD domain protein 1A XLP SAP MTCP1 LYP IMD5 EBVS DSHP RP5-1052M9.3 SH2 domain protein 1A, Duncan\1s disease (lymphoproliferative syndrome) SH2D1A
Names	
Usage	has to fail as it should be around 15kDa. [See for e.g., Czar (PNAS, Jun 2001; 98: 7449 - 7454) and Sayos (Blood, Jun 2001; 97: 3867 - 3874)]
Summary	
Accession ID	NP_002342.1
Immunogen	Peptide with sequence GIREDPDVCLKAP, from the C Terminus of the protein sequence according to NP_002342.1.
Peptide Sequence	GIREDPDVCLKAP
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	4068
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:1000.
Western Blot	Western blot: Approx 30-35kDa band seen in Jurkat cell lysate. Recommended for use at 2-4µg/ml.