

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

## **GOAT ANTI-SHB ANTIBODY**

**SKU:** EB05214

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

**Unit Size** 100 µg

Storage

Aliquot and store at -20°C. Minimize freezing and thawing. Instructions

Synonym /

SHB (Src homology 2 domain containing) adaptor protein B|bA3J10.2|RP11-3J10.8|RP11-3J10.3|Src homology 2

**Alias** domain containing adaptor protein BISHB **Names** 

**Accession** 

NP 003019.2

Peptide with sequence C-EHLSLLYPVAVRTL, from the C Terminus of the protein sequence according to **Immunogen** 

NP 003019.2.

**Peptide** 

C-EHLSLLYPVAVRTL Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

Method using the immunizing peptide.

Shipping

Refrigerated Instructions

**Predicted** 

Human, Mouse, Rat, Dog **Species** 

Reactive

Human **Species** 

Human

6461 **Gene ID** 

Mouse

230126

Gene ID

**Product** Grade

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

**ELISA** 

**Detection** Antibody detection limit dilution 1:6000.

Limit

Western Western blot: Approx 50kDa and 40kDa bands seen in Jurkat cell lysates. Recommended for use at 1-3µg/ml.

**Blot** Overnight incubation is recommended with this antibody.

Application

Pep-ELISA, WB **Type** 





Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600