Email: <a href="mailto:customerservice@vectorlabs.com">customerservice@vectorlabs.com</a>

Telephone: (650) 697-3600

## **GOAT ANTI-SHF ADAPTER PROTEIN ANTIBODY**

**SKU:** EB05160

×

## **SPECIFICATIONS**

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

 $\textbf{Unit Size} \qquad 100 \; \mu \text{g}$ 

Storage

Instructions

Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

**Alias** Shf adapter protein

Names

Accession Q7M4L6

\_. ..

Blocking EBP05160

Peptide

Immunogen
Peptide with sequence VAQSLSDAHHALTI-C, from the N Terminus of the protein sequence according to

Q7M4L6.

Peptide

Sequence VAQSLSDAHHALTI-C

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography

**Method** using the immunizing peptide.

Shipping

Instructions Re

Refrigerated

Predicted

Species Human

Reactive Species

Human

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite\_medium.png

Grade Western

Blot
Western blot: Approx 50kDa band seen in Human Brain extracts. Recommended for use at 0.5-2μg/ml.

Application

Pep-ELISA, WB

Type