

GOAT ANTI-SHF ADAPTER PROTEIN ANTIBODY

SKU: EB05160

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	Shf adapter protein
Accession ID	Q7M4L6
Blocking Peptide	EBP05160
Immunogen	Peptide with sequence VAQSLSDAHHALTI-C, from the N Terminus of the protein sequence according to Q7M4L6.
Peptide Sequence	VAQSLSDAHHALTI-C
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 50kDa band seen in Human Brain extracts. Recommended for use at 0.5-2µg/ml.
Application Type	Pep-ELISA, WB