

GOAT ANTI-SIX3B (ZEBRAFISH) ANTIBODY

SKU: EB07961



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	six3b sine oculis homeobox homolog 3b SIX3 cb347 six3.2 six6 sine oculis homeobox homolog 6
Accession ID	NP_571438.1
Blocking Peptide	EBP07961
Immunogen	Peptide with sequence C-NRLQHHGLGQSSLR, from the internal region of the protein sequence according to NP_571438.1 .
Product Comments	This antibody is not expected to cross-react with Six3a.
Peptide Sequence	C-NRLQHHGLGQSSLR
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Zebrafish
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:32000.
Western Blot	Not yet tested. At this stage we are dependent on researchers in the field for further characterization of this product. Therefore we cannot recommend an optimal concentration and the product is investigative grade. We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)