

# GOAT ANTI-SLAP2 ANTIBODY

**SKU:** EB05279



---

## SPECIFICATIONS

<b>Formulation</b>	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
<b>Unit Size</b>	100 µg
<b>Storage Instructions</b>	Aliquot and store at -20°C. Minimize freezing and thawing.
<b>Synonym / Alias Names</b>	Src-like adapter protein-2 src-like-adaptor 2 d 977B1.1 C20orf156 FLJ21992 SLAP-2 SLAP2 Src-like-adaptor 2 SLA2
<b>Accession ID</b>	NP_115590; NP_778252
<b>Blocking Peptide</b>	EBP05279
<b>Immunogen</b>	Peptide with sequence C-ISLNDEAVSLDDA, from the C Terminus of the protein sequence according to NP_115590; NP_778252.
<b>Peptide Sequence</b>	C-ISLNDEAVSLDDA
<b>Purification Method</b>	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
<b>Shipping Instructions</b>	Refrigerated
<b>Predicted Species</b>	Human
<b>Reactive Species</b>	Human
<b>Human Gene ID</b>	84174
<b>Product Grade</b>	<a href="https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png">https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png</a>
<b>Western Blot</b>	Western blot: Approx 35kDa band seen in MOLT-4 cell lysate. Recommended for use at 2-4µg/ml.