

GOAT ANTI-SLC16A4 ANTIBODY

SKU: EB05133



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	monocarboxylate transporter 4 MCT4 solute carrier family 16 (monocarboxylic acid transporters), member 4 SLC16A4
Accession ID	NP_004687
Blocking Peptide	EBP05133
Immunogen	Peptide with sequence C-FVPLAERWKNSLT, from the C Terminus of the protein sequence according to NP_004687.
Peptide Sequence	C-FVPLAERWKNSLT
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human
Reactive Species	Human
Human Gene ID	9122
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
Western Blot	Western blot: Approx 55-60kDa band seen in 293 cell lysate. Recommended for use at 0.5-2 µg/ml.