

GOAT ANTI-SLC9A1/ NHE-1 ANTIBODY

SKU: EB09311

250kDa

150kDa

100kDa

75kDa

50kDa

37kDa

25kDa

20kDa

15kDa

SPECIFICATIONS

Unit Size	100 µg
Synonym / Alias Names	solute carrier family 9, isoform A1 solute carrier family 9 (sodium/hydrogen exchanger), isoform 1 (antiporter, Na ⁺ /H ⁺ , amiloride sensitive) sodium/hydrogen exchanger 1 OTTHUMP00000004468 Na-Li countertransporter Na ⁺ /H ⁺ antiporter, amiloride-sensitive NHE1 FLJ42224 APNH solute carrier family 9 (sodium/hydrogen exchanger), member 1 (antiporter, Na ⁺ /H ⁺ , amiloride sensitive) SLC9A1
Accession ID	NP_003038.2
Blocking Peptide	EBP09311
Immunogen	Peptide with sequence C-RQKARQLEQKINNY , from the internal region of the protein sequence according to NP_003038.2.
Product Comments	The immunizing peptide represents part of the cytoplasmic domain
Peptide Sequence	C-RQKARQLEQKINNY
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Cow, Pig
Reactive Species	Human
Human Gene ID	6548
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:8000.
Western Blot	Approx 80kDa band observed in Human Kidney lysates (calculated MW of 90.8kDa according to NP_003038.2). Recommended concentration: 1-3µg/ml.
Application Type	Pep-ELISA, WB

GALLERY IMAGES

