


GOAT ANTI-SLIMMER / FHL1 ISOFORM B ANTIBODY

SKU: EB06508



250kDa
150kDa
100kDa
75kDa
50kDa
37kDa
25kDa
20kDa
15kDa
10kDa

SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	skeletal muscle LIM-protein1 Isl-1 and Mec-3 domains 1 four-and-a-half Lin11 LIM protein SLIMMER XMPMA SLIMMER, MGC111107 KYOT FLH1A FHL1B FHL1A RP11-535K18.1 bA535K18.1 (LIM protein SLIMMER) Four-and-a-half LIM domains 1 four and a half LIM domains 1 bA535K18.1 SLIM1 KYO-T FHL1
Accession ID	AAC72886.1
Blocking Peptide	EBP06508
Immunogen	Peptide with sequence CTVSRVSHPVSKARK, from the internal region of the protein sequence according to AAC72886.1.
Product Comments	This antibody is specific for the FHL-1 isoform B (SLIMMER) and is expected to recognise an epitope corresponding to aa 233-246 of human SLIMMER protein.
Peptide Sequence	CTVSRVSHPVSKARK
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Dog
Reactive Species	Human
Human Gene ID	2273
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/elite_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:128000.
Western Blot	Approx 30kDa band observed in Human Heart lysates is consistent with literature (PMID:10480922; calculated MW of 36.3kDa according to AAC72886). Recommended concentration: 0.01-0.1µg/ml.
Application Type	Pep-ELISA, WB

DOCUMENTS

- [Data Sheet](#)

GALLERY IMAGES

