

Email: customerservice@vectorlabs.com

Telephone: (650) 697-3600

GOAT ANTI-SLITRK1 ANTIBODY

SKU: EB07047

SPECIFICATIONS

Formulation Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.

Unit Size

Storage

Instructions Aliquot and store at -20°C. Minimize freezing and thawing.

Synonym /

Alias

SLITRK1|SLIT and NTRK-like family, member 1|HGNC:20297|KIAA0918|KIAA1910|LRRC12|RP11-395N17.1|leucine

rich repeat containing 12|slit and trk like 1 protein|slit and trk like gene 1|FL|54428

Names

Accession NP 443142.1

Blocking

EBP07047 **Peptide**

Peptide with sequence NADGAHRVYDCGSHS, from the C Terminus of the protein sequence according to Immunogen

NP 443142.1.

Peptide

NADGAHRVYDCGSHS Sequence

Purification Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using

Method the immunizing peptide.

114798

Shipping

Refrigerated Instructions

Predicted

Species

Human, Mouse, Rat, Dog, Pig, Cow

Human

Gene ID

Product

https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring medium.png

Grade **ELISA**

Detection

Antibody detection limit dilution 1:64000.

Limit

Preliminary experiments in human brain lysates gave no specific signal but low background (at antibody

Western concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been Blot

reported with other antibodies/lysates?

Application

Type

Pep-ELISA