

GOAT ANTI-SLITRK1 ANTIBODY

SKU: EB07047



SPECIFICATIONS

Formulation	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH7.3 with 0.5% bovine serum albumin.
Unit Size	100 µg
Storage Instructions	Aliquot and store at -20°C. Minimize freezing and thawing.
Synonym / Alias Names	SLITRK1 SLIT and NTRK-like family, member 1 HGNC:20297 KIAA0918 KIAA1910 LRRC12 RP11-395N17.1 leucine rich repeat containing 12 slit and trk like 1 protein slit and trk like gene 1 FLJ54428
Accession ID	NP_443142.1
Blocking Peptide	EBP07047
Immunogen	Peptide with sequence NADGAHRVYDCGSHS, from the C Terminus of the protein sequence according to NP_443142.1.
Peptide Sequence	NADGAHRVYDCGSHS
Purification Method	Purified from goat serum by ammonium sulphate precipitation followed by antigen affinity chromatography using the immunizing peptide.
Shipping Instructions	Refrigerated
Predicted Species	Human, Mouse, Rat, Dog, Pig, Cow
Human Gene ID	114798
Product Grade	https://prod-vector-labs-pimcore-assets.s3.us-east-1.amazonaws.com/assets/products/image/aspiring_medium.png
ELISA Detection Limit	Antibody detection limit dilution 1:64000.
Western Blot	Preliminary experiments in human brain lysates gave no specific signal but low background (at antibody concentration up to 1µg/ml). We would appreciate any feedback from people in the field - have any results been reported with other antibodies/lysates?
Application Type	Pep-ELISA

DOCUMENTS

- [Data Sheet](#)